

TECHNICAL CHART HI-INDEX 1.74 ASPHERIC

CYL (- ON -) - 0.00 ~ - 4.00					POWER		
DIA- (mm)	E/T (mm)	C/T (mm)	BACK BASE	FRONT BASE			
/					0.00		
					0.25		
					0.50		
					0.75		
					1.00		
					1.25		
					1.50		
					1.75		
					2.00		
Φ75	3.2	1.0	4.50	2.25	2.25		
	3.4			2.00	2.50		
	3.6			1.75	2.75		
	3.9			1.50	3.00		
	4.1	1.0	5.25	2.00	3.25		
	4.3			1.75	3.50		
	4.5			1.50	3.75		
	4.7			1.25	4.00		
	5.0	1.0	6.00	1.75	4.25		
	5.2			1.50	4.50		
	5.4			1.25	4.75		
	5.7			1.00	5.00		
	6.2	1.0	6.75	1.50	5.25		
	6.4			1.25	5.50		
	6.6			1.00	5.75		
	6.8			0.75	6.00		
	7.2	1.0	7.50	1.25	6.25		
	7.4			1.00	6.50		
	7.7			0.75	6.75		
	7.9			0.50	7.00		
8.2	1.0	8.25	1.00	7.25			
8.4			0.75	7.50			
8.6			0.50	7.75			
8.8			0.25	8.00			
Φ70 CYL- 0.00 ~ - 3.00	7.9	1.0	9.25	1.00	8.25		
	8.1			0.75	8.50		
	8.3			0.50	8.75		
	8.5			0.25	9.00		
	8.9			10.25	1.00	9.25	
	9.1				0.75	9.50	
	9.3				0.50	9.75	
9.5	0.25	10.00					
Φ70 CYL- 0.00 ~ - 2.00	10.2	1.0	11.00	0.50	10.50		
	10.7			0.50	11.00		
	11.3			0.50	11.50		
	11.9			0.50	12.00		
Φ65 CYL- 0.00	10.3	1.0	13.00	0.50	12.50		
	10.8			0.50	13.00		
	11.4			0.50	13.50		
	11.9			0.50	14.00		
	12.4			0.50	14.50		
12.9	0.50	15.00					

1. "MINUS" LENS
 SPH AXIS : -2.25 ~ -8.00
 CYL AXIS : -0.00 ~ -4.00
 & SPH AXIS: -8.25 ~ -10.00
 CYL AXIS: -0.00 ~ -3.00
 & SPH AXIS : -10.50 ~ -12.00
 CYL AXIS : -0.00 ~ -2.00
 & SPH AXIS : -12.50 ~ 15.00
 CYL AXIS : -0.00
2. "PLUS" LENS - NONE
3. TOLERANCE (+ 0.3/ - 0.1 mm)
 CENTER-E