

TECHNICAL CHART

HI-INDEX 1.74 ASPHERIC

CYL (- ON -) - 0.00 ~ - 3.00					POWER	
DIA- (mm)	E/T (mm)	C/T (mm)	BACK BASE	FRONT BASE		
<div style="position: absolute; top: 0; left: 0; width: 100%; height: 100%; border: 1px solid black; transform: rotate(45deg);"></div>					0.00	
					0.25	
					0.50	
					0.75	
					1.00	
					1.25	
					1.50	
					1.75	
					2.00	
Φ75	3.2	1.0	4.50	2.25	2.25	<ol style="list-style-type: none"> 1. "MINUS" LENS SPH AXIS : -2.25 ~ -10.00 CYL AXIS : -0.00 ~ -3.00 & SPH AXIS : -10.50 ~ -12.00 CYL AXIS : -0.00 ~ -2.00 & SPH AXIS : -12.50 ~ 15.00 CYL AXIS : -0.00 2. "PLUS" LENS - NONE 3. TOLERANCE (+ 0.3 / - 0.1 mm) CENTER-E
	3.4			2.00	2.50	
	3.6			1.75	2.75	
	3.9			1.50	3.00	
	4.1	1.0	5.25	2.00	3.25	
	4.3			1.75	3.50	
	4.5			1.50	3.75	
	4.7			1.25	4.00	
	5.0	1.0	6.00	1.75	4.25	
	5.2			1.50	4.50	
	5.4			1.25	4.75	
	5.7			1.00	5.00	
	6.2	1.0	6.75	1.50	5.25	
	6.4			1.25	5.50	
	6.6			1.00	5.75	
	6.8			0.75	6.00	
	7.2	1.0	7.50	1.25	6.25	
	7.4			1.00	6.50	
	7.7			0.75	6.75	
	7.9			0.50	7.00	
8.2	1.0	8.25	1.00	7.25		
8.4			0.75	7.50		
8.6			0.50	7.75		
8.8			0.25	8.00		
Φ70	7.9	1.0	9.25	1.00	8.25	
	8.1			0.75	8.50	
	8.3			0.50	8.75	
	8.5			0.25	9.00	
	8.9		10.25	1.00	9.25	
	9.1			0.75	9.50	
	9.3			0.50	9.75	
	9.5			0.25	10.00	
Φ70 CYL - 0.00 ~ - 2.00	10.2	1.0	11.00	0.50	10.50	
	10.7		11.50	0.50	11.00	
	11.3		12.00	0.50	11.50	
	11.9		12.50	0.50	12.00	
Φ65 CYL -0.00	10.3	1.0	13.00	0.50	12.50	
	10.8		13.50	0.50	13.00	
	11.4		14.00	0.50	13.50	
	11.9		14.50	0.50	14.00	
	12.4		15.00	0.50	14.50	
	12.9		15.50	0.50	15.00	