

TECHNICAL CHART

POLY ASP ULTRACLEAR EZ-EDGE AR

CYL (- ON -) - 0.00 ~ - 3.00					POWER	CYL (+ ON -) - 0.00 ~ - 3.00				
DIA- (mm)	E/T (mm)	C/T (mm)	BACK BASE	FRONT BASE		FRONT BASE	BACK BASE	C/T (mm)	E/T (mm)	DIA- (mm)
Φ72	2.0	2.0	3.85	3.83	0.00					
	2.0	1.7	4.07	3.83	0.25	3.88	3.65	2.0	1.8	Φ68
	2.3		4.30	3.83	0.50	3.88	3.43	2.2		
	2.6		4.52	3.83	0.75	3.88	3.20	2.3		
	2.9		4.75	3.83	1.00	3.88	2.97	2.5	1.6	
	3.2		4.98	3.83	1.25	3.88	2.75	2.6		
	3.5		1.4	5.20	3.83	1.50	3.88	2.53	2.8	
	3.4	3.61			1.75	3.88	2.33	2.6		
	3.6	3.38			2.00	3.88	2.10	2.9		
	3.9	3.16			2.25	4.20	2.22	2.9		
	4.2	2.93		2.50	4.20	2.00	3.1	1.0		
	4.5	3.16		2.75	4.20	1.78	3.4			
	4.8	2.93		3.00	4.20	1.56	3.6			
	5.0	2.70		3.25	4.20	1.33	3.9			
	5.2	1.3	5.65	2.48	3.50	5.28	2.20	4.0	-0.00 ~ -2.00	
	5.5			2.26	3.75	5.28	1.98	4.3		
	5.8			2.03	4.00	5.28	1.76	4.5		
	5.9			1.80	4.25	5.55		4.0		
	6.3	1.2	6.10	2.03	4.50	5.78	1.80	4.3	1.0	
	6.5			1.80	4.75	6.00		4.4		
	6.8			1.58	5.00	5.78	1.35	4.6		
	7.0			1.35	5.25	6.00		4.8		
	7.3			1.13	5.50	6.22		5.0		
	7.5			0.90	5.75	6.00	0.90	5.2		
7.7	0.68			6.00	6.22	5.4				
Φ72 CYL -0.00 ~ -2.00	8.3			1.2	6.56	0.90	6.25	1. "MINUS" LENS SPH: -0.00 ~ -6.00/-8.00 CYL: -0.00 ~ -3.00/-2.00 2. "PLUS" LENS SPH: +0.25 ~ +3.00/+6.00 CYL: -0.00 ~ -3.00/-2.00 3. TOLERANCE C/T: +/- 0.3mm Diameter: +1.00mm/-0.5mm		
	8.5	0.68	6.50							
	8.8	0.46	6.75							
	9.0	0.23	7.00							
	9.5	7.01	0.46		7.25					
	9.7		0.23		7.50					
	10.2	7.47	0.46		7.75					
	10.4		0.23		8.00					
-0.00 ~ -2.00										Φ60 CYL -0.00 ~ -2.00