

TECHNICAL CHART POLYCARBONATE EZ-TINT(GEN.2)

CYL (- ON -) - 0.00 ~ - 2.00					POWER	CYL (+ ON -) - 0.00 ~ - 2.00				
DIA- (mm)	E/T (mm)	C/T (mm)	BACK BASE	FRONT BASE		FRONT BASE	BACK BASE	C/T (mm)	E/T (mm)	DIA- (mm)
Φ70	2.20	2.20	4.45	4.45	0.00					
	2.30	2.00		4.20	0.25	2.25	2.17	1.73	1.50	Φ65
	2.50	1.90		3.95	0.50	2.40		1.95		
	2.60	1.80		3.70	0.75	2.50		2.18		
	2.90			3.45	1.00	2.80		2.11		
	3.00	1.60		3.20	1.25	3.10		2.34		
	3.20			2.95	1.50	3.25		2.56		
	3.40	1.40	5.95	4.20	1.75	3.50		2.78		
	3.60			3.95	2.00	3.75	3.02			
	3.90			3.70	2.25	3.50	3.30			
	4.20			3.45	2.50	3.75	3.47			
	4.40			3.20	2.75	4.00	3.70			
	4.70			2.95	3.00	4.25	3.92			
	5.00			1.30	7.45	4.20	3.25	4.50	1.71	
	5.30	3.95	3.50			4.25	4.38			
	5.60	3.70	3.75			4.75	4.61			
	5.90	3.45	4.00			5.00	4.94			
	6.10	3.20	4.25							
	6.40	2.95	4.50							
	6.80	7.95	3.20		4.75					
	7.00		2.95		5.00					
	7.30		2.70		5.25					
	7.60		2.45		5.50					
	7.80		2.20	5.75						
	8.10	1.95	6.00							