

TECHNICAL CHART POLYCARBONATE EZ-TINT(GEN.2)

CYL (- ON -) - 0.00 ~ - 2.00					POWER	CYL (+ ON -) - 0.00 ~ - 2.00				
DIA- (mm)	E/T (mm)	C/T (mm)	BACK BASE	FRONT BASE		FRONT BASE	BACK BASE	C/T (mm)	E/T (mm)	DIA- (mm)
Spherical Φ70	2.20	2.20	4.45	4.45	0.00					
	2.30	2.00		4.20	0.25	2.25	2.17			
	2.50	1.90		3.95	0.50	2.40		1.95		
	2.60	1.80		3.70	0.75	2.50		2.18		
	2.90			3.45	1.00	2.80		2.11		
	3.00	1.60		3.20	1.25	3.10		2.34		
	3.20			2.95	1.50	3.25		2.56		
	3.40	1.40	5.95	4.20	1.75	3.50		2.78	1.20	
	3.60			3.95	2.00	3.75	3.02			
	3.90			3.70	2.25	3.50	3.30			
	4.20			3.45	2.50	3.75	3.47			
	4.40			3.20	2.75	4.00	3.70			
	4.70			2.95	3.00	4.25	3.92			
	5.00			1.30	7.45	4.20	3.25	4.50		1.71
	5.30	3.95	3.50			4.25	4.38			
	5.60	3.70	3.75			4.75	4.61			
	5.90	3.45	4.00			5.00	4.94			
	6.10	3.20	4.25							
	6.40	2.95	4.50							
	6.80	3.20	4.75							
7.00	2.95	5.00								
7.30	2.70	5.25								
7.60	2.45	5.50								
7.80	2.20	5.75								
8.10	1.95	6.00								