

TECHNICAL CHART POLYCARBONATE SUNSMART AR

CYL (- ON -) - 0.00 ~ - 2.00					POWER	CYL (+ ON -) - 0.00 ~ - 2.00				
DIA- (mm)	E/T (mm)	C/T (mm)	BACK BASE(1.586)	FRONT BASE(1.586)		FRONT BASE(1.586)	BACK BASE(1.586)	C/T (mm)	E/T (mm)	DIA- (mm)
Spherical Φ70	2.20	2.20	4.46	4.45	0.00					
	2.30	2.00		4.20	0.25					
	2.50	1.90		3.95	0.50	2.89	2.40	1.70	1.50	
	2.60	1.80		3.70	0.75	3.14		1.91		
	2.90			3.45	1.00	3.39		2.14		
	3.00	1.60		3.20	1.25	3.63		2.07		
	3.20			2.95	1.50	3.88		2.30		
	3.40	1.40	5.94	4.20	1.75	4.11		2.55		
	3.60			3.95	2.00	4.36		2.76		
	3.90			3.70	2.25	4.11	2.98			
	4.20			3.45	2.50	4.36	3.22			
	4.40			3.20	2.75	4.59	3.44			
	4.70			2.95	3.00	4.83	3.56			
	5.00			1.89	7.44	4.20	3.25	5.07	3.82	
	5.30	3.95	3.50			5.30	4.06			
	5.60	3.70	3.75			5.54	4.25			
	5.90	3.45	4.00			5.88	4.49			
	6.10	3.20	4.25							
	6.40	2.95	4.50							
	6.80	3.20	4.75							
7.00	2.95	5.00								
7.30	2.70	5.25								
7.60	2.45	5.50								
7.80	2.20	5.75								
8.10	1.95	6.00								
Aspheric Φ65										