

TECHNICAL CHART

CR-39 SPHERICAL SEMI H/C

70mm Diameter Lens Blank

BASE	NOMINAL	TRUE CURVE (Nd=1.530)		RADIUS	SAG AT 48 (48mm)		THICKNESS (mm)		COMMENTS
	FRONT	FRONT	BACK	FRONT	FRONT	BACK	CENTER	EDGE	
2.50	2.50	2.13	6.12	248.860	1.16	3.39	9.6	14.5	
4.00	4.00	3.57	6.12	148.590	1.95	3.39	7.8	12.0	
5.00	5.00	5.59	5.00	94.800	3.09	2.75	10.6	9.8	
6.00	6.00	6.05	6.12	87.540	3.35	3.39	7.3	7.4	