

**TECHNICAL CHART**  
**POLYCARBONATE PROGRESSIVE ASP SEMI H/C(EZ-VIEW MINI)**

**75mm Diameter Lens Blank**

BASE	TRUE CURVE (Nd=1.530)		RADIUS	SAG AT 50 (50mm)		THICKNESS (mm)		ADDS IN 0.25 STEPS	RX Range Diopter
	FRONT	BACK	FRONT	FRONT	BACK	CENTER	EDGE		
<b>2.75</b>	2.93	6.25	180.887	1.74	3.77	7.2	13.2	1.00 to 3.00	-8.00 ~ -2.75
<b>4.75</b>	4.85	6.25	109.278	2.90	3.77	7.4	9.2	1.00 to 3.00	-2.50 ~ +1.75
<b>6.75</b>	6.76	6.25	78.402	4.09	3.77	8.0	8.8	1.00 to 3.00	+2.00 ~ +3.50
<b>7.75</b>	7.71	6.25	68.742	4.71	3.77	10.5	6.8	1.00 to 3.00	+3.75 ~ +8.00