

**TECHNICAL CHART**  
**POLYCARBONATE PROGRESSIVE ASP SEMI H/C(EZ-VIEW STD.)**

**75mm Diameter Lens Blank**

BASE	TRUE CURVE (Nd=1.530)		RADIUS	SAG AT 50 (50mm)		THICKNESS (mm)		ADDS IN 0.25 STEPS	RX Range  Diopter
	FRONT	BACK	FRONT	FRONT	BACK	CENTER	EDGE		
<b>3.50</b>	3.31	6.25	160.121	1.96	3.77	8.0	13.2	1.00 to 3.00	-7.00 ~ 0.00
<b>5.50</b>	5.20	6.25	101.923	3.11	3.77	8.5	9.2	1.00 to 3.00	0.00 ~ +4.00
<b>7.50</b>	7.18	6.25	73.816	4.36	3.77	9.6	6.8	1.00 to 3.00	+3.00 ~ +7.00