

**NEW**  
**SEMI-FINISHED**  
**PHOTOCHROMIC-GRAY**

## POLYCARBONATE SPHERICAL SFSV PHOTOCROMIC-GRAY

### SOMO Spherical Semi-Finished Single Vision Polycarbonate Hardcoat Photochromic-gray

Activation Speed highly comparable to Transitions-VI

Deactivation speed is *FASTER* than Transitions-VI


Significant cost saving vs. Transitions-VI

Premium Polycarbonate Resin

“Best in Class” Hard Coat Chemistry for photochromic polycarbonate

SOMO quality validated with COLTS independent testing

Spherical Base Curve specifications with Free Form tolerance




### Specifications

- ◆ Spherical Design
- ◆ 1.586 Refractive Index
- ◆ 1.20 Specific Gravity
- ◆ 31 Abbe Value
- ◆ 99.5% Visible Light Transmittance
- ◆ 100% UV Protection at 380nm
- ◆ Super Scratch Resistance Hardcoat
- ◆ 75 mm Blank Size

**BASE | 2.00 | 4.00 | 6.00 | 8.00**



**IT'S ALL ABOUT YOU**



**PHOTOCROMIC-BROWN  
COMING SOON  
SUMMER 2011**