

SEMI-FINISHED PHOTOCHROMIC-GRAY

SOMO Spherical Semi-Finished Single Vision Polycarbonate Hardcoat Photochromic-gray

Activation Speed highly comparable to Transitions-VI

Deactivation speed is **FASTER** than Transitions-VI


Significant cost saving vs. Transitions-VI

Premium Polycarbonate Resin

“Best in Class” Hard Coat Chemistry for photochromic polycarbonate

SOMO quality validated with COLTS independent testing

Spherical Base Curve specifications with Free Form tolerance



Specifications

- ◆ Spherical Design
- ◆ 1.586 Refractive Index
- ◆ 1.20 Specific Gravity
- ◆ 31 Abbe Value
- ◆ 99.5% Visible Light Transmittance
- ◆ 100% UV Protection at 380nm
- ◆ Super Scratch Resistance Hardcoat
- ◆ 75 mm Blank Size

BASE | 2.00 | 4.00 | 6.00 | 8.00



IT'S ALL ABOUT YOU