

TECHNICAL CHART

HI-INDEX 1.60 SPHERICAL

CYL (- ON -) - 0.00 ~ - 4.00					POWER	CYL (+ ON -) - 0.00 ~ - 2.00					
DIA- (mm)	E/T (mm)	C/T (mm)	BACK BASE	FRONT BASE		FRONT BASE	BACK BASE	C/T (mm)	E/T (mm)	DIA- (mm)	
Φ75	2.1	1.9	5.50	5.50	0.00						
	2.3	1.8		5.25	0.25	4.75	4.50	2.2	2.1	Φ70	
	2.4	1.6		5.00	0.50	5.00		2.1	1.7		
	2.6			4.75	0.75	5.25		2.3			
	2.9			4.50	1.00	5.50		2.3	1.4		
	3.1	1.4		4.25	1.25	5.75		2.5	1.5		
	3.3			4.00	1.50	6.00		2.8			
	3.5			3.75	1.75	6.25		3.2			
	3.6			3.50	2.00	6.50		3.5			
	4.0	1.1	6.50	4.25	2.25	5.75		3.50	3.2		1.4
	4.3			4.00	2.50	6.00	3.5				
	4.6			3.75	2.75	6.25	3.7				
	4.9			3.50	3.00	6.50	4.0				
	5.2			3.25	3.25	6.75	4.3		1.5		
	5.6			3.00	3.50	7.00	4.5				
	5.9			2.75	3.75	7.25	4.7				
	6.2			2.50	4.00	7.50	4.9				
	6.7	1.1	7.50	3.25	4.25	5.75	1.50	4.9	1.4		
	6.9			3.00	4.50	6.00		5.1			
	7.1			2.75	4.75	6.25		5.3			
	7.5			2.50	5.00	6.50		5.6			
	7.8			2.25	5.25	6.75		5.9	1.5		
	8.1			2.00	5.50	7.00		6.1			
	8.3			1.75	5.75	7.25		6.3			
	8.6			1.50	6.00	7.50		6.6			
	Φ70 CYL -0.00 ~ -2.00	8.2	1.1	8.50	2.25	6.25	1. "MINUS" LENS SPH AXIS : -0.00 ~ -6.00/-8.00 CYL AXIS : -0.00 ~ -4.00/-2.00 2. "PLUS" LENS SPH AXIS : +0.25 ~ +6.00 CYL AXIS : -0.00 ~ -2.00 3. TOLERANCE (± 0.3 mm)				
		8.5			2.00	6.50					
		8.7			1.75	6.75					
9.0		1.50			7.00						
9.2		1.25			7.25						
9.5		1.00			7.50						
9.7		0.75			7.75						
10.0		0.50			8.00						