

1.74 & Poly- tintable®

• For the Active, Outdoors Type



Not only is sunlight damaging to the skin, it also impairs the eye and can even cause cataracts and corneal burns. Enjoy the freedom of the outdoors without the risk of exposing your eyes to damaging light.

• For Drivers on the Road



Lights from vehicles on the opposite lane, reflected light and sunlight are all factors that block our vision. SOMO's tintable lenses are the ideal choice in protecting your eyes from these dangerous glares.

The World's Very First



Pending patent registration, SOMO Optical's proprietary technology brings the world a 1.74 tintable lens.

Freedom in Color Selection



Can tint in a variety of colors to match the customer's own swatch and taste.

Enhanced Impact Resistance



Protect your eyes from potential risks with FDA approved impact resistance.

Long Lasting Visual Clarity



Experience long lasting, comfortable vision with visual acuity, lens durability and improved light transmittance by combining SOMO's state-of-the-art product quality and coating technology.

1.74 Spherical SF Single Vision - Hardcoated

Base (nominal)	Diameter (mm)	Thickness (mm)		RX Range
		Center	Edge	
1.00	75	7.5	14.0	-16.00D to -6.25D
2.00		5.3	11.5	-6.00D to -2.75D
3.00		6.0	10.5	-2.50D to +1.00D
4.00		6.0	9.0	+1.25D to +2.00D
5.00		8.5	9.5	+2.25D to +3.50D
6.00		8.7	8.4	+3.75D to +5.00D
7.00		9.0	7.2	+5.25D to +6.50D
8.00		11.0	7.7	+6.75D to +8.00D



1.59 Polycarbonate Spherical SF Single Vision - Hardcoated*

Base (nominal)	Diameter (mm)	Thickness (mm)		RX Range
		Center	Edge	
5.00	80	7.5	13.0	-10.00D to -2.25D
6.00		6.5	15.1	-10.00D to -6.25D
7.00		8.0	14.8	-10.00D to -6.25D
8.00		9.0	14.0	-10.00D to -6.25D

*Extra thick blanks

