

**SOMO** [ **ULTRACLEAR AR**  
with ITO ]

**ITO ANTI-STATIC LAYER**  
> reduces dust particles  
> easiest to clean

**OLEOPHOBIC/ HYDROPHOBIC TOPCOAT**  
on ULTRACLEAR AR with ITO  
> grease-resistant  
> easier to clean

**HYDROPHOBIC TOPCOAT**  
on Premium Clear AR with ITO  
> water-resistant  
> easy to clean

**SUPREME MULTI ANTI-REFLECTIVE LAYERS**  
> reduces reflections from computer screens and televisions  
> reduces eye fatigue  
> reduces distorting glare  
> effectively eliminates reflections on the lens surface


**INDEX MATCHED THERMAL HARD COAT LAYER**  
> toughest scratch resistance available

**ULTRA BONDING LAYER**  
> extends the life of the lens


**HIGH IMPACT CUSHION COATING**  
> provides spectacular durability

**LENS SUBSTRATE**  
> 1.74, 1.67, 1.60, 1.59 polycarbonate, and 1.50 CR-39 plastic

**THE COLTS PERFORMANCE SEAL**  
> validates that SOMO lenses will deliver consistent product performance and quality that you can depend on




**Top-rated SOMO [ ULTRACLEAR AR with ITO ] cost significantly less than the leading AR brands.**




**[ ULTRACLEAR AR with ITO ] Available in CR39, 1.60, 1.67, 1.74 and Polycarbonate**

The UltraClear high performance finished lens series has three top coats that provide excellent cleanability by repelling water, oils and with the ITO anti static layer, dust and dirt.

All UltraClear lenses have the COLTS Seal of Certifiable Quality and have earned "best in class" performance ratings for the COLTS Real Life Simulation Wearers Test.



Somo Quartz Coatings      Other Coatings



Other Coatings      Ultraclear ITO      Other Coatings      Ultraclear ITO

Somo Optical's Quartz AR coating process makes lenses scratch and smudge resistant, prevents moisture buildup and repels contaminants making it easier to keep lenses clean. Includes ITO anti-static layer.

**ENHANCED AR TECHNOLOGY AT THE BEST VALUE.**



www.somoptical.com