

TECHNICAL CHART

1.74 ASPHERIC UV420+ BLUE BLOCK ULTRACLEAR AR

CYL (- ON -) - 0.00 ~ - 2.50					POWER	
DIA- (mm)	E/T (mm)	C/T (mm)	BACK BASE	FRONT BASE		
<div style="position: relative; width: 100%; height: 100%; border: 1px solid black;"> </div>					0.00	
					0.25	
					0.50	
					0.75	
					1.00	
					1.25	
					1.50	
					1.75	
					2.00	
					2.25	
					2.50	
					2.75	
					3.00	
					3.25	
					3.50	
3.75						
4.00						
Φ75	5.0	1.0	6.00	1.75	4.25	1. "MINUS" LENS SPH AXIS : -4.00 ~ -10.00 CYL AXIS : -0.00 ~ -2.50 & SPH AXIS : -10.50 ~ -11.00 CYL AXIS : -0.00 ~ -2.00 2. "PLUS" LENS - NONE 3. TOLERANCE (+ 0.3 / - 0.1 mm) CENTER-E
	5.2			1.50	4.50	
	5.4			1.25	4.75	
	5.7			1.00	5.00	
	6.2	1.0	6.75	1.50	5.25	
	6.4			1.25	5.50	
	6.6			1.00	5.75	
	6.8			0.75	6.00	
	7.2	1.0	7.50	1.25	6.25	
	7.4			1.00	6.50	
	7.7			0.75	6.75	
	7.9			0.50	7.00	
	8.2	1.0	8.25	1.00	7.25	
	8.4			0.75	7.50	
	8.6			0.50	7.75	
8.8	0.25			8.00		
Φ70	7.9	1.0	9.25	1.00	8.25	
	8.1			0.75	8.50	
	8.3			0.50	8.75	
	8.5			0.25	9.00	
	8.9		10.25	1.00	9.25	
	9.1			0.75	9.50	
	9.3			0.50	9.75	
	9.5			0.25	10.00	
Φ70 CYL 0.00 ~ 2.00	10.2	1.0	11.00	0.50	10.50	
	10.7		11.50	0.50	11.00	