

# TECHNICAL CHART

## 1.67 ASP SUNSMART PHOTO-GRAY PREMIUM AR

CYL ( - ON - ) - 0.00 ~ - 2.00					POWER	CYL ( + ON - ) - 0.00 ~ - 2.00					
DIA- (mm)	E/T (mm)	C/T (mm)	BACK BASE	FRONT BASE		FRONT BASE	BACK BASE	C/T (mm)	E/T (mm)	DIA- (mm)	
$\Phi 75$	3.40	1.10	4.50	2.50	2.00	1. "MINUS" LENS  SPH AXIS : -2.00 ~ -6.00  CYL AXIS : -0.00 ~ -2.00  2. TOLERANCE ( $\pm 0.3$ mm )  CEN'YL -0.00 ~ -1.00					
	3.60			5.00	2.75						2.25
	3.80				2.50						2.50
	4.00		2.25		2.75						
	4.20		2.00		3.00						
	4.60		5.50		2.25						3.25
	4.70			2.00	3.50						
	4.90			1.75	3.75						
	5.20			1.50	4.00						
	5.60		6.00	1.75	4.25						
	5.80			1.50	4.50						
	6.00			1.25	4.75						
	6.20			1.00	5.00						
	6.50		6.50	1.25	5.25						
	6.80			1.00	5.50						
	7.10			0.75	5.75						
	7.30			0.50	6.00						