

TECHNICAL CHART

POLYCARBONATE UV420+ BLUE BLOCK ULTRACLEAR AR

CYL (- ON -) - 0.00 ~ - 2.00					POWER	CYL (+ ON -) - 0.00 ~ - 2.00						
DIA- (mm)	E/T (mm)	C/T (mm)	BACK BASE	FRONT BASE		FRONT BASE	BACK BASE	C/T (mm)	E/T (mm)	DIA- (mm)		
$\Phi 70$	2.20	2.20	4.03	4.03	0.00	/						
	2.30	2.00		3.80	0.25	2.17	2.39	1.70	1.50	$\Phi 65$		
	2.50	1.90		3.57	0.50			2.62			1.91	
	2.60	1.80		3.35	0.75			2.84			2.14	
	2.90			3.12	1.00			3.06			2.07	
	3.00	1.60		2.89	1.25			3.28	2.30			
	3.20			2.67	1.50			3.51	2.55			
	3.40	1.40		5.37	3.80			1.75	3.73		2.76	
	3.60		3.57		2.00			3.94	2.98			
	3.90		3.35		2.25	3.77	3.22					
	4.20		3.12		2.50	3.93	3.44					
	4.40		2.89		2.75	4.15	3.56					
	4.70		2.67		3.00	4.37	3.82					
	5.00	1.30	6.73	3.80	3.25	4.58	1.71	4.06				
	5.30			3.57	3.50	4.80		4.25				
	5.60			3.35	3.75	5.01		4.49				
	5.90			3.12	4.00	5.32		5.04				
	6.10		7.19	7.19	2.89	4.25		1.71	1. Plus lenses are referred to edge thickness at spherical power edge. 2. Center thickness tolerance shall be less than +/- 0.3mm 3. Diameter tolerance shall be less than +1.0mm / -0.5mm			
	6.40				2.67	4.50						
	6.80				2.89	4.75						
7.00	2.67				5.00							
7.30	2.44	5.25										
7.60	2.22	5.50										
7.80	1.99	5.75										
8.10	1.76	6.00										