



Poly SunColors® Sun-curve Pre-Tinted FSV

Design Features

- First poly Sun-Curve design in FSV category 3 sunglasses.
- 80mm large size & cover 60% of popular base 6 sunglasses frame.
- 100% UV400 block sunglasses.
- Provide FSV backside premium AR coating
- High scratch resistance hard coat properties than classical cooked tinted lens .
- Tinted Gray is natural and unadulterated perception
- Visual light transmission : 15%.
- Available: FSV +2.00 to -4.00 with cylinder up to -2.00



Advantages :

- Provide consumers want prescription with their fashion sunglasses frame .
- Eliminate hassle of polycarbonate lens tinting process.
- Good solution to optical retailer to offer less one-day prescription special treatment sunglasses to their clients.
- Available replace one lens as it have color consistency to match the old perfectly.