

TECHNICAL CHART

CR-39 UV420+ BLUE BLOCK ULTRACLEAR AR

CYL (- ON -) - 0.00 ~ - 2.00					POWER	CYL (+ ON -) - 0.00 ~ - 2.00				
DIA- (mm)	E/T (mm)	C/T (mm)	BACK BASE	FRONT BASE		FRONT BASE	BACK BASE	C/T (mm)	E/T (mm)	DIA- (mm)
Φ70	2.0	2.0	5.00	5.00	0.00					
	2.3			4.75	0.25	5.25				
	2.6			4.50	0.50	5.50				
	2.8	1.9		4.25	0.75	5.75	5.00	2.2	1.9	
	3.2			4.00	1.00	6.00		2.5		
	3.4			3.75	1.25	5.75		2.7		
	3.7	1.8		3.50	1.50	6.00	4.50	2.8	1.6	
	4.0			3.25	1.75	6.25		3.0		
	4.3			3.00	2.00	6.50		3.1	1.4	
	4.6			2.75	2.25	6.25		3.2		
	4.9	1.7	2.50	2.50	6.50	4.00	3.5	1.2		
	5.2		2.25	2.75	6.75		3.7			
	5.5		2.00	3.00	7.00		4.0			
	5.8		1.75	3.25	6.25		4.3			
	6.1		1.50	3.50	6.50		4.7			
	6.7		2.75	3.75	6.75		4.9			
	6.9	6.50	2.50	4.00	7.00	3.00	5.1			
	7.2		2.25	4.25			5.4			
	7.5		2.00	4.50			5.8			
	8.0		2.25	4.75						
8.3	1.6	7.00	2.00	5.00						
8.6			1.75	5.25						
8.9			1.50	5.50						
9.2			1.25	5.75						
9.5			1.00	6.00						

1. TOLERANCE OF FRONT/BACK CURVE: ±0.25
2. C/T & E/T TOLERANCE: ± 0.2
3. POWER TOLERANCE: S -0.00 to S -4.00 0.09D
 S -4.25 to S -6.00 0.12D
 S -0.00/C -2.00 to S -4.00/C -2.00 0.09D
 S -4.25/C -2.00 to S -6.00/C -2.00 0.12D