

TECHNICAL CHART FSV POLY UV420+ PHOTOCROMIC- GY ULTRACLEAR AR
--

CYL (- ON -) - 0.00 ~ - 2.00					POWER	CYL (+ ON -) - 0.00 ~ - 2.00					
DIA- (mm)	E/T (mm)	C/T (mm)	BACK BASE	FRONT BASE		FRONT BASE	BACK BASE	C/T (mm)	E/T (mm)	DIA- (mm)	
Spherical Φ70	2.20	2.20	4.03	4.03	0.00						
	2.30	2.00		3.80	0.25	2.39	2.17	1.70	1.50	Aspheric Φ65	
	2.50	1.90		3.57	0.50	2.62		1.91			
	2.60	1.80		3.35	0.75	2.84		2.14			
	2.90			3.12	1.00	3.06		2.07	1.30		
	3.00	1.60		2.89	1.25	3.28		2.30			
	3.20			2.67	1.50	3.51		2.55			
	3.40	1.40	5.37	3.80	1.75	3.73		2.76			
	3.60			3.57	2.00	3.94		2.98			
	3.90			3.35	2.25	3.77	1.71	3.22			
	4.20			3.12	2.50	3.93		3.44			
	4.40			2.89	2.75	4.15		3.56			
	4.70			2.67	3.00	4.37		3.82			
	5.00	1.30	6.73	3.80	3.25	4.58		4.06			
	5.30			3.57	3.50	4.80		4.25			
	5.60			3.35	3.75	5.01		4.49			
	5.90			3.12	4.00	5.32		5.04			
	6.10			2.89	4.25	1. "MINUS" LENS SPH: -0.00 ~ -6.00 CYL: -0.00 ~ -2.00 2. "PLUS" LENS SPH: +0.25 ~ +4.00 CYL: -0.00 ~ -2.00 3. TOLERANCE C/T: +/- 0.3mm Diameter: +1.00mm/-0.5mm					
	6.40			2.67	4.50						
	6.80		7.19	2.89	4.75						
	7.00			2.67	5.00						
	7.30			2.44	5.25						
	7.60			2.22	5.50						
	7.80			1.99	5.75						
	8.10			1.76	6.00						